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EU Type examination carried out by : APAVE SUDEUROPE SAS -
 CS60193 - 13322 MARSEILLE CEDEX 16 – France

Manufacturing inspection carried out by : QUINTIN Certifications
 (n 2927) 825 Route de Romans - 38160 St-Antoine-l'Abbaye - France

Reference	Class
2120L	FFP2 NR
2130	FFP3 NR
2130V	FFP3 NR

Filtering half masks to protect against particles EN149 : 2001 + A1 : 2009

In compliance with PPE regulation 2016/425

Category III PPE

EU declarations of conformity available at : www.inspire-protection.fr/declarations-de-conformite/


WARNING – This notice must be read and understood before using equipment.

The purpose of the single-use respiratory mask that you have acquired is to protect you against solid and liquid particles, as well as non-toxic fibers. A mask certified EN149 : 2001 + A1 : 2009 must be classified in one of three classes FFP1 - FFP2 or FFP3. Only the marking on the mask indicates with certainty the class of protection of the mask.

1. You must not use this respiratory mask without first defining the environment where you will be using it so as to validate whether it specifically corresponds to your requirements. Should these instructions not be complied with, you may be putting your health, or even your life, at risk. Do not hesitate to call upon a professional (e.g. a safety officer or occupational physician) to help you define the protection you require.
2. Not respecting the instructions herein may compromise your health. Please respect them carefully.
3. As this mask does not have an oxygen supply, you must use it in a sufficiently ventilated area. If the oxygen level is less than 17%, you will rapidly experience dizziness and breathing difficulties. It is imperative to return to a well-ventilated area in order to avoid more acute symptoms.
4. Never use this mask if you do not know the precise nature of the toxic materials you may come into contact with. In particular, do not use it if you are in an area containing explosive products or atmospheres. This mask does not protect against gases and fumes.
5. The mask's protection is less effective when worn over a beard. Indeed, facial hair makes it difficult to ensure that the mask is airtight, particularly around the edges. Please be extra careful when adjusting it on your face to make sure that it is impossible to breathe in air which has not been filtered through the mask. Please refer to the below instructions before use.
6. 2130 FFP3 NR mask fits better on medium faces. Please consult the below size table for face dimensions recommendations. The right fit of the mask can vary from one user to another, so please carefully adjust to guarantee proper seal on the face before use. Please refer to the below instructions before use.
7. Never attempt to modify the mask which has been designed and certified for your safety (e.g. do not change the elastic straps, do not write on the mask, do not change the clip etc.)
8. This mask is for single use (the marking NR on the mask : Not Re-usable). After a maximum 8 hours of continuous use, destroy it. No maintenance operation will enable it to recover its protective capacities after use. Washing the mask is pointless and may even be dangerous for you. Do not put back on your face a mask which has been used (even for less than 8 hours) and taken off. Remove it also if you have difficulty breathing (breathlessness, headache), regardless of how long you have been wearing it.
9. Keep your masks in a place where the maximum humidity is 70% at a temperature between -25 C and +40 C. so as to maintain its characteristics. After x years following the date of manufacture (see manufacture date and use-by date on the packaging as well as on the PPE (label EXP YYYY/MM on the mask), please consider the mask to be out-of-date, it can no longer be safely used, and destroy it.
10. In any case and in spite of all the care taken during the manufacture of this respiratory protection mask, it is the responsibility of the user alone (and in no case that of the manufacturer) to ensure the adequation of the mask with the characteristics of the environment where it is used. The manufacturer shall, in no case, guarantee any specific length of life for the mask, due to their lack of knowledge of the environment where the mask has been stored after leaving the factory and of the environment where it will be used (temperature, pressure, hygrometry, particle concentration to be filtered, user's morphology, physical effort etc.)
11. This mask shall protect you correctly only if you have carefully followed these recommendations as well as the instructions for use outlined hereafter:

For FFP3 NR ref 2130	A+B+C
Maximum face dimensions (mm)	380

Facial dimensions

HOW TO USE THE MASK

Check that the nose clip is present and that it does not pierce the filter. Check that there are two elastic straps properly sealed to the mask. Hold the mask and open the two folds. Shape the nose clip using your thumb and index fingers. With the other hand, take hold of the two elastic straps, stretch and pass them behind your head. The lower strap needs to be positioned under your ears, and the upper strap needs to be placed on top of your cranium, above your ears. Take hold of the second strap and pass it over the top of your head and place it above your ears. Adjust the mask on your face, placing the nose clip against your nose. Using your thumb and index finger, pinch the nose clip so that it is snug over your nose and does not let any air come between the mask and your face. Take the time to fit the edges of the mask to every contour of your face so that no unfiltered air can pass under it. If the mask is not sufficiently airtight you can shorten the elastic straps by tying a knot in them to make the mask airtight. The care you take at this stage will determine the real effectiveness of the mask and therefore the user's health. Please pay the utmost attention to putting the mask on correctly. Do not hesitate to adjust the elastic straps if necessary. These masks have been designed to allow the wearing of glasses : we recommend that the glasses be placed over the mask at the level of the nose clip. Finally, place your hands on the edges of the mask and breathe out. If you feel no air, then your mask has been fitted properly and you can start using it. If you feel air coming out the edges, re-adjust using the method described above one time only. If you do not manage to make the mask airtight, remove it and dispose of it.

WHEN TO DISPOSE OF YOUR MASK?

After any use (end of shift, removal of mask, demonstration of how to fit mask etc.)
 If you notice any defects in the mask (e.g. the nose clip is not attached or has pierced the fabric, defective or missing elastic straps, mask edges not welded together properly, non-conform marking)
 If you do not manage to make the mask airtight after a maximum of two attempts
 If you experience dizziness, breathlessness, headaches or any other more acute symptoms while wearing the mask
 If a defect occurs while you are wearing the mask (e.g. the elastic gives)

Meaning of the pictograms on packaging

REF Product reference

LOT Manufacture batch N 

Manufacturing details DD/MM/YY HH:MM are specified on each mask